



# **3-D Rotating Mixer 40 rpm**

## **LX911RM**

## 3-D Rotating Mixer 40 rpm LX911RM

3-D Rotating Mixer 40 rpm LX911RM is a mixer with horizon, vertical and angled rotation with a speed variance of 10 rpm ~ 40 rpm. Adjustable tube holder makes the product either horizontally or vertically or between the angles of horizon-vertical to rotate. Equipped with 2 adjustable rotating frames to meet different hybrid modes in one mixture.

### Features

- Combination of adjustable rotating frame
- Variance of rotation modes - mild type to somersault
- Reciprocating oscillation function from 20° to 99°
- Provision with three types of rotating frame

### Applications

Used in molecular hybridization, blood banking and various samples requiring mixing

### Specifications

Model No	LX911RM
Speed range ( fixed )	10 rpm ~ 40 rpm
Reciprocating swing angle range	20 ° ~ 99 °
Time	1 s ~ 999 min
Voltage	AC 220 V / AC 110 V ; 50 / 60 Hz
Power	4 W
Fuse	250 V ; 1 A ; ϕ5 x 20
Dimension ( W x D x H )	240 x 145 x 210 mm
Weight	1.2 kgs

### Optional Accessories

Accessory No	Standard configuration	Description
RM01	1 pair ( 2 pieces )	5 centrifuge tubes 10 ml / 15 ml ( ϕ15 mm ~ ϕ16 mm ) + 6 centrifuge tubes 5 ml / 7 ml ( ϕ12 mm ~ ϕ13 mm )
RM02	1 pair ( 2 pieces )	3 centrifuge tubes 50 ml + 18 centrifuge tubes 0.5 ml
RM03	1 pair ( 2 pieces )	18 centrifuge tubes 1.5 ml / 2.0 ml