



36 L Muffle Furnace

LX669MF

www.labdex.com info@labdex.com

36 L Muffle Furnace LX669MF

36 L Muffle Furnace LX669MF is a microprocessor controlled chamber with stable PID temperature control system. Equipped with a rapid heating speed of 0 to 30 °C / min, the chamber adopts vacuum forming ceramic fiber with no power drops off with high temperature. Temperature uniformity with excellent protection of heating elements (SiC rod) is maintained. One press start button for pre-set programmes with alerting signal for overheating of chamber.

Features

- High temperature heating alloy wire for heating elements
- Rapid heating speed 0 to 30 °C / min
- PID controlled temperature controller with a precision of 1 °C
- One button press start for pre-programmed functions
- Durable inner lining of light insulating bricks
- Double walled housing for high stability
- 1700 °C MoSi₂ heating rod
- Alert signal for overheating of chamber
- Heat curve can be preset with the programmable function
- 2 programmable groups with 8 segments in each group

Application

Used for fusing glass, soldering, brazing, creating enamel coatings, ceramics, research centers, medical laboratories to determine the non-volatile and non-combustible proportion of the sample and for highly sophisticated metallurgical applications

Specifications

LX669MF
36 L
1400 °C
1350 °C
SiC
± 1 °C
0 to 30 °C / min
2
8
Alumina Ceramic Fiber
S style
380 V
10 kW
300 x 400 x 300 mm

Standard accessories

Accessory No	Accessory
1	Crucible tong
2	High temperature gloves

Optional Accessories

Accessory No	Accessory
1	Paperless recorder
2	Alumina crucible
3	Stainless steel exhaust chamber