



# Digital Rocking Shaker

LX138DRS

www.labdex.com info@labdex.com

## **Digital Rocking Shaker LX138DRS**

Digital Rocking Shaker LX138DRS is a natural rubber and stainless steel platform. It is a non-deforming, easy to clean and prevention unit. Equipped with real-time display of speed and diminishing time. With a speed range of 10 rpm  $\sim$  80 rpm and 10  $^{\circ}$  swing motion this unit is ideal for laboratory use.

#### **Features**

- Stainless steel and natural rubber platform which is alkali-resistant, non-deforming and easy to clean preventing spillage
- Speed and diminishing time display
- Articulated rotating bearing for prolonged use

# **Applications**

Used for low foaming, agitating and DNA extractions, staining and de-staining procedures

### **Specifications**

Model No	LX138DRS
Speed range	10 rpm ~ 80 rpm
Swing cycloid	25 mm ( up and down ) / 10 $^{\rm 0}$
Time range	1 min ~ 99 h 59 min or continuous
Maximum load	5.0 kgs ( double - deck platform )
Platform dimension	228 x 280 mm
Fuse	250 V, 1A, \( 5 \times 20
Voltage	AC 100 ~ 230 V; 50 / 60 Hz
Power	50 W
Dimension (W x D x H)	284 x 264 x 248 mm
Weight	6.0 kgs